# WELLSACLEAN B

## BOHRERBAD

REF WCC605

### Produkttyp

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REF WCC605

### Product description

WELLSACLEAN B BOHRERBAD is a ready to use drill disinfection solution. It contains potassium hydroxide, alcohols and amines as well as corrosion inhibitors for material protection. The extensive spectrum of efficacy and the very good material compatibility offer a safety level for the applicant during the reprocessing and application of rotating precision instruments. WELLSACLEAN B BOHRERBAD is free of aldehydes and phenols.

### Areas of application

WELLSACLEAN B BOHRERBAD is applied for the manual disinfection of rotating invasive medical devices (e.g. drills, mills, diamond grinders and other precision instruments) in an immersion bath.

### Suitable for:

- Stainless steel drills
- Diamond drills
- Mills
- Diamond grinders

Not suitable for aluminium parts, elastic polishers, ceramic grinding instruments and root canal instruments.

Do not use for squirting or spraying. Do not use for disinfection of skin and hands.

### Material compatibility

<b>Metals</b>	stainless steel
<b>Plastics</b>	PE, PSU, PVC
<b>Rubber</b>	silicone, EPDM, natural rubber
<b>Ceramic</b>	porcelain
<b>Advice</b>	Upon compliance to the intended use and the indicated contact time, no material damage is to be expected with the above mentioned materials. <b>Avoid contact with aluminium and linoleum.</b> Test material compatibility for polycarbonate, polyacrylate, brass, and copper before use.

### Application

Put on protective clothing, protective gloves, and safety glasses prior to disinfection! Do not allow debris to dry on the instruments.

### Disinfection procedure:

1. Pour undiluted drill disinfection solution into the drill bat. Immerse instruments. Instruments must be completely covered by the solution.
2. After 5 minutes of contact time: Remove instruments from the drill bath and rinse them thoroughly with water (at least drinking water quality, approx. 20 °C) and dry them.
3. In case of persistent debris: Brush off debris under running water, repeat the disinfection procedure with 5 min contact time in fresh drill disinfection solution.
4. Afterwards, check instruments for cleanliness and integrity. Comply with the instrument manufacturer's recommendations for reprocessing as well as with national guidelines for processing of medical devices.

### Shelf life and expiry

Shelf life: Renew daily and in case of visible contamination. The expiry date of the product is indicated on the label. Do not use the product after the expiry date.

### Special Advice

**Report serious incidents related to this product to the manufacturer and the respective authority.**

**Do not mix the product with other chemicals. Prior to the first application of the product or after the application of another product, the containers have to be cleaned and rinsed thoroughly.**

**Rinse the instruments with demineralized water if they are supposed to be sterilized.**

**Application only by qualified personnel.**

**Do not pour disinfection solution that has been removed back into the original container.**

# WELLSACLEAN B

## BOHRERBAD

REF WCC605

### Labelling according to Regulation (EC) No 1272/2008 (CLP)

For information on labelling according to Regulation (EC) No 1272/2008 (CLP) see safety datasheet and label.

### Spectrum of efficacy

Efficacy under dirty conditions (undiluted, at room temperature)	Standard	Contact time
<b>BACTERIA</b>		
S. aureus, E. hirae, P. aeruginosa	EN 13727	1 min
	EN 14561	1 min
S. aureus, E. hirae, E. coli, P. mirabilis, P. aeruginosa	Standard methods VAH	5 min
<b>YEASTS</b>		
C. albicans	EN 13624	1 min
	EN 14562	1 min
C. albicans	Standard methods VAH	5 min
<b>FUNGI</b>		
A. brasiliensis	Standard methods VAH	5 min
<b>MYCOBACTERIA/TB</b>		
M. terrae, M. avium	EN 14348	2 min
	EN 14563	5 min
M. terrae, M. avium	Standard methods VAH	5 min
<b>VIRUSES</b>		
MVA modified Vacciniavirus Ankara	EN 14476	1 min
Polyomavirus SV40		5 min
Poliovirus		
Adenovirus		
MNV murine Norovirus		

### Expert opinions and additional information

THE EFFICACY OF THE WELLSACLEAN B BOHRERBAD has been assessed by independent laboratories.

The listings, reports, and notes in this product information refer to the product PLIWA Roto-Pren, which is identical in formulation.

For further details concerning expert opinions, information on products, and answers to specific questions contact us by e-mail (info@wellsamed.com).

The manufacturer has a manufacturing authorisation according as per §13 German medicines act as well as a GMP- and an EN 13485 certificate.

### Identification, certificates

CE 0123 MD IIb

The PLIWA Roto-Pren, which has the same formulation as the WELLSACLEAN B BOHRERBAD, is VAH-listed.

### Product composition

100 g ready to use solution contain:

15.00 g Isopropanol, 2.97 g potassium hydroxide, 2.10 g N,N-Bis(3-aminopropyl)-dodecylamine (> 90 %), 1.25 g sodium nitrite, 1.12 g borax decahydrate

### Physico-chemical data

appearance: clear, colourless liquid  
 pH-value: > 13  
 density (20 °C): 1.010 – 1.030 g/cm<sup>3</sup>

### Storage advice

Store between 5 and 25 °C. Keep the container tightly closed in a well-ventilated place. Store product in an upright position. Keep away from heat, hot surfaces, sparks, open flames, and other ignition sources. No smoking. Protect from sunlight. Keep only in original container.

See safety datasheet for further advice.

 PLIWA Hygiene GmbH,  
 Pliwastraße 2, 34323 Malsfeld, Germany

Distributed by:  
 wellsamed GmbH  
 Brahestraße 16, 04347 Leipzig, Germany